




SLC6A7 Antibody

Product Code	CSB-PA04209A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99884
Immunogen	Recombinant Human Sodium-dependent proline transporter protein (557-636AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sodium-dependent proline transporter (Solute carrier family 6 member 7), SLC6A7, PROT
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SLC6A7

Image


Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA04209A0Rb at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.