



SLC6A7 Antibody, HRP conjugated

Product Code	CSB-PA04209B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99884
Immunogen	Recombinant Human Sodium-dependent proline transporter protein (557-636AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sodium-dependent proline transporter (Solute carrier family 6 member 7), SLC6A7, PROT
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SLC6A7