




ANGPTL3 Antibody

Product Code	CSB-PA05759A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5C1
Immunogen	Recombinant Human Angiotensin-related protein 3 protein (17-268AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Angiotensin-related protein 3 (Angiotensin-5) (ANG-5) (Angiotensin-like protein 3) [Cleaved into: ANGPTL3(17-221); ANGPTL3(17-224)], ANGPTL3, ANGPT5
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	ANGPTL3

Image



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA05759A0Rb at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.