



CASP1 Antibody, HRP conjugated

Product Code	CSB-PA05829B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29466
Immunogen	Recombinant Human Caspase-1 protein (120-269AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Caspase-1 (CASP-1) (EC 3.4.22.36) (Interleukin-1 beta convertase) (IL-1BC) (Interleukin-1 beta-converting enzyme) (ICE) (IL-1 beta-converting enzyme) (p45) [Cleaved into: Caspase-1 subunit p20; Caspase-1 subunit p10], CASP1, IL1BC IL1BCE
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CASP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.