





MMP2 Antibody

Product Code	CSB-PA06879A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08253
Immunogen	Recombinant Human 72 kDa type IV collagenase protein (445-660AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	72 kDa type IV collagenase (EC 3.4.24.24) (72 kDa gelatinase) (Gelatinase A) (Matrix metalloproteinase-2) (MMP-2) (TBE-1) [Cleaved into: PEX], MMP2, CLG4A
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	MMP2


Image



Immunohistochemistry of paraffin-embedded human endometrial cancer using CSB-PA06879A0Rb at dilution of 1:100



Immunohistochemistry of paraffin-embedded human prostate cancer using CSB-PA06879A0Rb at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA06879A0Rb at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.