



CD14 Antibody, HRP conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA07085B0Rb |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P08571 |
| Immunogen | Recombinant Human Monocyte differentiation antigen CD14 protein (20-345AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Monocyte differentiation antigen CD14 (Myeloid cell-specific leucine-rich glycoprotein) (CD antigen CD14) [Cleaved into: Monocyte differentiation antigen CD14, urinary form; Monocyte differentiation antigen CD14, membrane-bound form], CD14 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Immunology |
| Target Names | CD14 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |