



CD14 Antibody, Biotin conjugated

Product Code	CSB-PA07085D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08571
Immunogen	Recombinant Human Monocyte differentiation antigen CD14 protein (20-345AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Monocyte differentiation antigen CD14 (Myeloid cell-specific leucine-rich glycoprotein) (CD antigen CD14) [Cleaved into: Monocyte differentiation antigen CD14, urinary form; Monocyte differentiation antigen CD14, membrane-bound form], CD14
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CD14
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.