




Thpo Antibody

Product Code	CSB-PA07649A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49745
Immunogen	Recombinant Rat Thrombopoietin protein (22-194AA)
Raised In	Rabbit
Species Reactivity	Rat, Mouse, Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thrombopoietin, Thpo
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Cardiovascular
Target Names	Thpo

Image



Western blot
 All lanes: Thpo antibody at 4μg/ml
 Lane 1: Mouse liver tissue
 Lane 2: 293T whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 35 kDa
 Observed band size: 35 kDa