



Thpo Antibody, HRP conjugated

Product Code	CSB-PA07649B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49745
Immunogen	Recombinant Rat Thrombopoietin protein (22-194AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thrombopoietin, Thpo
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Cardiovascular
Target Names	Thpo
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.