



CXCL3 Antibody, Biotin conjugated

Product Code	CSB-PA08999D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19876
Immunogen	Recombinant Human C-X-C motif chemokine 3 protein (35-107AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	C-X-C motif chemokine 3 (GRO-gamma(1-73)) (Growth-regulated protein gamma) (GRO-gamma) (Macrophage inflammatory protein 2-beta) (MIP2-beta) [Cleaved into: GRO-gamma(5-73)], CXCL3, GRO3 GROG SCYB3
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CXCL3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.