



MYL2 Antibody, Biotin conjugated

Product Code	CSB-PA10059D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10916
Immunogen	Recombinant Human Myosin regulatory light chain 2, ventricular/cardiac muscle isoform protein (8-164AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myosin regulatory light chain 2, ventricular/cardiac muscle isoform (MLC-2) (MLC-2v) (Cardiac myosin light chain 2) (Myosin light chain 2, slow skeletal/ventricular muscle isoform) (MLC-2s/v) (Ventricular myosin light chain 2), MYL2, MLC2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MYL2