



NRG4 Antibody, Biotin conjugated

Product Code	CSB-PA10569D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WWG1
Immunogen	Recombinant Human Pro-neuregulin-4, membrane-bound isoform protein (1-62AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pro-neuregulin-4, membrane-bound isoform (Pro-NRG4) [Cleaved into: Neuregulin-4 (NRG-4)], NRG4
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NRG4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.