



GKN1 Antibody, HRP conjugated

Product Code	CSB-PA10589B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NS71
Immunogen	Recombinant Human Gastrokine-1 protein (49-199AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gastrokine-1 (18 kDa antrum mucosa protein) (AMP-18) (Protein CA11), GKN1, AMP18 CA11
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	GKN1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.