





HSPA12B Antibody

Product Code	CSB-PA10989A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96MM6
Immunogen	Recombinant Human Heat shock 70 kDa protein 12B protein (201-500AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Heat shock 70 kDa protein 12B, HSPA12B, C20orf60
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	HSPA12B

Image



Immunohistochemistry of paraffin-embedded human endometrial cancer using CSB-PA10989A0Rb at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA10989A0Rb at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.