



HSPA12B Antibody, Biotin conjugated

Product Code	CSB-PA10989D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96MM6
Immunogen	Recombinant Human Heat shock 70 kDa protein 12B protein (201-500AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Heat shock 70 kDa protein 12B, HSPA12B, C20orf60
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	HSPA12B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.