



LRG1 Antibody, Biotin conjugated

Product Code	CSB-PA11059D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02750
Immunogen	Recombinant Human Leucine-rich alpha-2-glycoprotein protein (36-346AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leucine-rich alpha-2-glycoprotein (LRG), LRG1, LRG
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	LRG1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.