



TG Antibody, FITC conjugated

Product Code	CSB-PA12325C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01266
Immunogen	Recombinant Human Thyroglobulin protein (2500-2767AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thyroglobulin (Tg), TG
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.