



ACTA2 Antibody, HRP conjugated

Product Code	CSB-PA12659B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62736
Immunogen	Recombinant Human Actin, aortic smooth muscle protein (3-265AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Actin, aortic smooth muscle (Alpha-actin-2) (Cell growth-inhibiting gene 46 protein) [Cleaved into: Actin, aortic smooth muscle, intermediate form], ACTA2, ACTSA ACTVS
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ACTA2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.