



RNASE3 Antibody, HRP conjugated

Product Code	CSB-PA12789B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12724
Immunogen	Recombinant Human Eosinophil cationic protein (28-160AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Eosinophil cationic protein (ECP) (EC 3.1.27.-) (Ribonuclease 3) (RNase 3), RNASE3, ECP RNS3
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RNASE3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.