



ACTG2 Antibody, FITC conjugated

Product Code	CSB-PA12799C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P63267
Immunogen	Recombinant Human Actin, gamma-enteric smooth muscle protein (181-376AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Actin, gamma-enteric smooth muscle (Alpha-actin-3) (Gamma-2-actin) (Smooth muscle gamma-actin) [Cleaved into: Actin, gamma-enteric smooth muscle, intermediate form], ACTG2, ACTA3 ACTL3 ACTSG
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ACTG2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.