



ATG5 Antibody, HRP conjugated

Product Code	CSB-PA15499B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H1Y0
Immunogen	Recombinant Human Autophagy protein 5 protein (5-272AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Autophagy protein 5 (APG5-like) (Apoptosis-specific protein), ATG5, APG5L ASP
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ATG5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.