



PRKCH Antibody, HRP conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA15879B0Rb |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P24723 |
| Immunogen | Recombinant Human Protein kinase C eta type protein (1-678AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Protein kinase C eta type (EC 2.7.11.13) (PKC-L) (nPKC-eta), PRKCH, PKCL PRKCL |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Signal Transduction |
| Target Names | PRKCH |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |