



NEU1 Antibody, Biotin conjugated

Product Code	CSB-PA17249D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99519
Immunogen	Recombinant Human Sialidase-1 protein (201-415AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sialidase-1 (EC 3.2.1.18) (Acetylneuraminyl hydrolase) (G9 sialidase) (Lysosomal sialidase) (N-acetyl-alpha-neuraminidase 1), NEU1, NANH
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NEU1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.