



MT1B Antibody, HRP conjugated

Product Code	CSB-PA17419B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07438
Immunogen	Recombinant Human Metallothionein-1B protein (1-61AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Metallothionein-1B (MT-1B) (Metallothionein-IB) (MT-IB), MT1B, MT1Q
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MT1B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.