



MT4 Antibody, FITC conjugated

Product Code	CSB-PA17459C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P47944
Immunogen	Recombinant Human Metallothionein-4 protein (1-60AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Metallothionein-4 (MT-4) (Metallothionein-IV) (MT-IV), MT4
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MT4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.