



choD Antibody, HRP conjugated

Product Code	CSB-PA18185B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P9WMV9
Immunogen	Recombinant Mycobacterium tuberculosis Cholesterol oxidase protein (1-578AA)
Raised In	Rabbit
Species Reactivity	Mycobacterium tuberculosis
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cholesterol oxidase (EC 1.1.3.6) (Cholesterol isomerase) (EC 5.3.3.1), choD
Immunogen Species	Mycobacterium tuberculosis (strain ATCC 25618 / H37Rv)
Research Area	Others
Target Names	choD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.