



MT2731 Antibody, Biotin conjugated

Product Code	CSB-PA300323LD01MVZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P9WJ10
Immunogen	Recombinant Mycobacterium tuberculosis Antitoxin MT2731 protein (1-81AA)
Raised In	Rabbit
Species Reactivity	Mycobacterium tuberculosis
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Antitoxin MT2731
Immunogen Species	Mycobacterium tuberculosis (strain CDC 1551 / Oshkosh)
Research Area	Others
Target Names	MT2731
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.