



cagA Antibody

Product Code	CSB-PA301307LA01HUS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P80200
Immunogen	Recombinant Helicobacter pylori Cytotoxicity-associated immunodominant antigen protein (918-1147AA)
Raised In	Rabbit
Species Reactivity	Helicobacter pylori
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytotoxicity-associated immunodominant antigen (120 kDa protein), cagA, cai
Immunogen Species	Helicobacter pylori (Campylobacter pylori)
Research Area	Others
Target Names	cagA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.