



cagA Antibody, HRP conjugated

Product Code	CSB-PA301307LB01HUS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P80200
Immunogen	Recombinant Helicobacter pylori Cytotoxicity-associated immunodominant antigen protein (918-1147AA)
Raised In	Rabbit
Species Reactivity	Helicobacter pylori
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytotoxicity-associated immunodominant antigen (120 kDa protein), cagA, cai
Immunogen Species	Helicobacter pylori (Campylobacter pylori)
Research Area	Others
Target Names	cagA