



III Antibody

Product Code	CSB-PA303116LA01ECY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P69168
Immunogen	Recombinant Enterobacteria phage M13 Attachment protein G3P protein (19-424AA)
Raised In	Rabbit
Species Reactivity	Enterobacteria phage M13
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Attachment protein G3P (Gene 3 protein) (G3P) (Minor coat protein), III
Immunogen Species	Enterobacteria phage M13 (Bacteriophage M13)
Research Area	Others
Target Names	III
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.