



Ribosome-inactivating protein PD-L1/PD-L2 Antibody

Product Code	CSB-PA306294LA01PJJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P84853
Immunogen	Recombinant Phytolacca dioica Ribosome-inactivating protein PD-L1/PD-L2 protein (1-261AA)
Raised In	Rabbit
Species Reactivity	Phytolacca dioica
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribosome-inactivating protein PD-L1/PD-L2 (EC 3.2.2.22) (rRNA N-glycosidase PD-L1/PD-L2)
Immunogen Species	Phytolacca dioica
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.