



Ribosome-inactivating protein cucurmosin Antibody, Biotin conjugated

Product Code	CSB-PA307485LD01CUI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P84531
Immunogen	Recombinant Cucurbita moschata Ribosome-inactivating protein cucurmosin protein (1-244AA)
Raised In	Rabbit
Species Reactivity	Cucurbita moschata
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribosome-inactivating protein cucurmosin (EC 3.2.2.22) (rRNA N-glycosidase)
Immunogen Species	Cucurbita moschata
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.