



Ribosome-inactivating protein lychnin Antibody, Biotin conjugated

Product Code	CSB-PA307648HD01LKN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P85101
Immunogen	Recombinant Silene chalcedonica Ribosome-inactivating protein lychnin protein (1-234AA)
Raised In	Rabbit
Species Reactivity	Silene chalcedonica
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribosome-inactivating protein lychnin (EC 3.2.2.22)
Immunogen Species	Silene chalcedonica
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.