



Ribosome-inactivating protein charybdin Antibody, FITC conjugated

Product Code	CSB-PA308125HC01DEZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P84786
Immunogen	Recombinant Drimia maritima Ribosome-inactivating protein charybdin protein (1-257AA)
Raised In	Rabbit
Species Reactivity	Drimia maritima
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribosome-inactivating protein charybdin (EC 3.2.2.22) (rRNA N-glycosidase)
Immunogen Species	Drimia maritima
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.