



Scytovirin Antibody, Biotin conjugated

Product Code	CSB-PA308458LD01SEH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P86041
Immunogen	Recombinant Scytonema varium Scytovirin protein (1-95AA)
Raised In	Rabbit
Species Reactivity	Scytonema varium
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Scytovirin (SVN)
Immunogen Species	Scytonema varium
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.