



# Pollen allergen Dac g 3 Antibody, Biotin conjugated

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-PA309006HD01DAC  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                    |
| <b>Uniprot No.</b>         | P93124   |
| <b>Immunogen</b>           | Recombinant Dactylis glomerata Pollen allergen Dac g 3 protein (1-96AA)          |
| <b>Raised In</b>           | Rabbit   |
| <b>Species Reactivity</b>  | Dactylis glomeRata   |
| <b>Tested Applications</b> | ELISA  |
| <b>Form</b>                | Liquid   |
| <b>Conjugate</b>           | Biotin   |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| <b>Purification Method</b> | >95%, Protein G purified   |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Alias</b>               | Pollen allergen Dac g 3 (Allergen Dac g III) (allergen Dac g 3) (Fragment)       |
| <b>Immunogen Species</b>   | Dactylis glomerata   |
| <b>Research Area</b>       | Others   |