



Pollen allergen Dac g 3 Antibody, Biotin conjugated

Product Code	CSB-PA309006HD01DAC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P93124
Immunogen	Recombinant Dactylis glomerata Pollen allergen Dac g 3 protein (1-96AA)
Raised In	Rabbit
Species Reactivity	Dactylis glomeRata
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pollen allergen Dac g 3 (Allergen Dac g III) (allergen Dac g 3) (Fragment)
Immunogen Species	Dactylis glomerata
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.