



GP Antibody, HRP conjugated

Product Code	CSB-PA310843LB01ZAA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P87671
Immunogen	Recombinant Zaire ebolavirus Envelope glycoprotein (502-637AA)
Raised In	Rabbit
Species Reactivity	Zaire ebolavirus
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Envelope glycoprotein (GP1,2) (GP) [Cleaved into: GP1; GP2; GP2-delta], GP
Immunogen Species	Zaire ebolavirus (strain Eckron-76) (ZEBOV) (Zaire Ebola virus)
Research Area	Others
Target Names	GP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.