



p46 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA313811HA01MWK |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P0C0J7 |
| Immunogen | Recombinant Mycoplasma hyopneumoniae 46 kDa surface antigen protein (28-416AA) |
| Raised In | Rabbit |
| Species Reactivity | Mycoplasma hyopneumoniae |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | 46 kDa surface antigen (p46), p46 |
| Immunogen Species | Mycoplasma hyopneumoniae (strain 232) |
| Research Area | Others |
| Target Names | p46 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |