



p46 Antibody, HRP conjugated

Product Code	CSB-PA313811HB01MWK
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C0J7
Immunogen	Recombinant Mycoplasma hyopneumoniae 46 kDa surface antigen protein (28-416AA)
Raised In	Rabbit
Species Reactivity	Mycoplasma hyopneumoniae
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	46 kDa surface antigen (p46), p46
Immunogen Species	Mycoplasma hyopneumoniae (strain 232)
Research Area	Others
Target Names	p46
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.