



GP Antibody, Biotin conjugated

Product Code	CSB-PA314099LD01ZAT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C773
Immunogen	Recombinant Zaire ebolavirus Super small secreted glycoprotein (33-298AA)
Raised In	Rabbit
Species Reactivity	Zaire ebolavirus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Super small secreted glycoprotein (SsGP), GP
Immunogen Species	Zaire ebolavirus (strain Kikwit-95) (ZEBOV) (Zaire Ebola virus)
Research Area	Others
Target Names	GP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.