



trxC Antibody, Biotin conjugated

Product Code	CSB-PA314858HD01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0AGG4
Immunogen	Recombinant Escherichia coli Thioredoxin 2 protein (1-139AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thioredoxin 2 (Trx-2) (EC 1.8.1.8) (Protein-disulfide reductase), trxC, yfiG
Immunogen Species	Escherichia coli
Research Area	Others
Target Names	trxC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.