



PMP2 Antibody, HRP conjugated

Product Code	CSB-PA315130LB01HO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C6G6
Immunogen	Recombinant Horse Myelin P2 protein (2-132AA)
Raised In	Rabbit
Species Reactivity	Horse
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myelin P2 protein, PMP2
Immunogen Species	Equus caballus (Horse)
Research Area	Others
Target Names	PMP2