



PMP2 Antibody, Biotin conjugated

Product Code	CSB-PA315130LD01HO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C6G6
Immunogen	Recombinant Horse Myelin P2 protein (2-132AA)
Raised In	Rabbit
Species Reactivity	Horse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myelin P2 protein, PMP2
Immunogen Species	Equus caballus (Horse)
Research Area	Others
Target Names	PMP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.