



lasA Antibody

Product Code	CSB-PA318588LA01EZ X
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14789
Immunogen	Recombinant Pseudomonas aeruginosa Protease LasA protein (237-418AA)
Raised In	Rabbit
Species Reactivity	Pseudomonas aeruginosa
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protease LasA (EC 3.4.24.-) (Staphylolytic protease), lasA
Immunogen Species	Pseudomonas aeruginosa (strain ATCC 15692 / DSM 22644 / CIP 104116 / JCM 14847 / LMG 12228 / 1C / PRS 101 / PAO1)
Research Area	Others
Target Names	lasA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.