



Pka-C1 Antibody, FITC conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA318954LC01DLU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P12370 |
| Immunogen | Recombinant Drosophila melanogaster cAMP-dependent protein kinase catalytic subunit 1 protein (2-353AA) |
| Raised In | Rabbit |
| Species Reactivity | Drosophila melanogaster |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | cAMP-dependent protein kinase catalytic subunit 1 (PKA C) (Protein kinase DC0), Pka-C1, CdkA DC0 |
| Immunogen Species | Drosophila melanogaster |
| Research Area | Others |
| Target Names | Pka-C1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |