



Pka-C1 Antibody, Biotin conjugated

Product Code	CSB-PA318954LD01DLU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12370
Immunogen	Recombinant Drosophila melanogaster cAMP-dependent protein kinase catalytic subunit 1 protein (2-353AA)
Raised In	Rabbit
Species Reactivity	Drosophila melanogaster
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	cAMP-dependent protein kinase catalytic subunit 1 (PKA C) (Protein kinase DC0), Pka-C1, CdkA DC0
Immunogen Species	Drosophila melanogaster
Research Area	Others
Target Names	Pka-C1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.