



AP24 Antibody, HRP conjugated

Product Code	CSB-PA320333LB01NHE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14170
Immunogen	Recombinant Nicotiana tabacum Osmotin protein (22-246AA)
Raised In	Rabbit
Species Reactivity	Nicotiana tabacum
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Osmotin, AP24
Immunogen Species	Nicotiana tabacum (Common tobacco)
Research Area	Others
Target Names	AP24
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.