



vpx Antibody

Product Code	CSB-PA320847HA01SHA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11266
Immunogen	Recombinant Simian immunodeficiency virus Protein Vpx protein (1-112AA)
Raised In	Rabbit
Species Reactivity	Simian immunodeficiency virus
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein Vpx (Viral protein X) (X ORF protein), vpx
Immunogen Species	Simian immunodeficiency virus (isolate K78) (SIV-mac) (Simian immunodeficiency virus rhesus monkey)
Research Area	Others
Target Names	vpx
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.