



Major pepsin inhibitor 3 Antibody, Biotin conjugated

Product Code	CSB-PA322527LD01DOS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19400
Immunogen	Recombinant Ascaris suum Major pepsin inhibitor 3 protein (21-169AA)
Raised In	Rabbit
Species Reactivity	Ascaris suum
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major pepsin inhibitor 3 (PI-3)
Immunogen Species	Ascaris suum (Pig roundworm) (Ascaris lumbricoides)
Research Area	Others